


Rotavapor® R-300

Technical data sheet

The Rotavapor® R-300 meets the highest expectations in convenience and versatility for rotary evaporation. Its modular design allows for easy extension of the R-300 to a fully integrated system where a central interface regulates each component.





Overview of Rotavapor® R-300 stand-alone and systems configurations

	Stand-alone configura- tion	System configuration
Rotavapor® R-300 / R-300 Hand Lift**	1	1
Heating Bath B-301 / B-305**	*	1
Interface I-300 / I-300 Pro (incl. VacuBox)**	*	1
Vacuum Pump V-300	-	1
Recirculating Chiller F-305 / F-308**	-	*

^{*} Optional, according to order code

Rotavapor® R-300

The Rotavapor® R-300 represents high reliability and convenience in laboratory rotary evaporation applications. The Rotavapor® R-300 is available with either an electronic or manual lift, and can optionally include the Heating Bath B-301 or B-305 and, a glass assembly of your choice and the Interface I-300 or I-300 Pro.

All accessories and connections are included for out of the box usage.

^{**} According to order code

Order Code

Choose the configuration according to your needs:

1 1 R 3 0 0

Lift System

- 1 Manual lift
- 2 Electronic lift

Heating Bath

- 0 Without
- 1 B-301: 20 95 °C (water), max. flask size: 1 L
- 5 B-305: 20 220 °C, max. flask size: 5 L

Standard Joint -

- 0 Without (w/o evaporating flask, w/o vapor duct)
- 1 SJ 29/32 (incl. evaporating flask, vapor duct)
- 2 SJ 24/40 (incl. evaporating flask, vapor duct)

Glass Assembly

- N Without (w/o receiving flask) 1
- V Vertical (incl. receiving flask)
- A Diagonal (incl. receiving flask)
- C Cold trap (incl. receiving flask)
- S Vertical reflux (incl. receiving flask)

Protective Coating

- 0 Without 1
- 1 Safety coating (P+G) ²

Interface

- 0 Without
- 1 Without, with Woulff bottle
- 2 I-300, with Woulff bottle (incl. VacuBox)
- 3 I-300 Pro, with Woulff bottle (incl. VacuBox)
- 4 I-300, with valve unit (incl. VacuBox)
- 5 I-300 Pro, with valve unit (incl. VacuBox)

Voltage -

- 1 220 240 V
- 2 100 120 V

¹ If glass assembly E, CR, HP or BY is required, for Glass Assembly select "N", for Standard Joint select either "1" or "2" and for Protective Coating select "0" (vapor duct and evaporating flask are included). Glass assembly E, CR, HP or BY need to be ordered separately.

² "Safety coating (P+G)" applies to condenser and receiving flask. PLASTIC+GLAS coating is sensitive to higher vapor temperatures.


Scope of delivery of Rotavapor® R-300 stand-alone configuration

All configurations are supplied ready to use.

Components	Stand-alone configura- tion
Rotavapor® R-300 / R-300 Hand Lift**	1
Heating Bath B-301 / B-305**	*
Glass assembly**	*
Vapor duct (with Combi-Clip)	**
Vacuum seal VS26	1
Evaporating flask (1 liter)	**
Receiving flask (1 liter)	**
Flyer "List of solvents"	1
Power cable (1 for heating bath, 1 for Rotavapor®)	**
Set of required cooling tubings	**
Interface I-300 / I-300 Pro**	*
VacuBox	*
Communications cable. BUCHI COM 6 p	*

^{*} Optional, according to order code

^{**} According to order code



Rotavapor® R-300 system configuration

The rotary evaporator R-300 system meets the highest expectations in convenience and versatility. It is bundled with the Vacuum Pump V-300, controlled by the Interface I-300 / I-300 Pro and can include the Recirculating Chiller F-305 / F-308. The interface is the central control unit for all process parameters. Rotation speed, heating-, cooling-, vapor-temperature and pressure are perfectly synchronized to optimize the distillation process while increasing efficiency and convenience. All the Rotavapor® R-300 system configurations are P+G coated. All accessories and connections are included for out of the box usage.

Choose the configuration according to your needs:

1 1 S R 3 0 0

Lift System

- 1 Manual lift
- 2 Electronic lift

Heating Bath

- 1 B-301: 20 95 °C (water), max. flask size: 1 L
- 5 B-305: 20 220 °C, max. flask size: 5 L

Standard Joint

- 1 SJ 29/32 (incl. evaporating flask, vapor duct)
- 2 SJ 24/40 (incl. evaporating flask, vapor duct)

Glass Assembly

- V Vertical (condenser and receiving flask are P+G coated)¹
- C Cold trap (condenser and receiving flask are P+G-LT coated)12

Interface

- S I-300, with Woulff bottle (incl. VacuBox)
- P I-300 Pro, with Woulff bottle (incl. VacuBox)

Recirculating Chiller

- 0 Without
- 1 F-305: cooling capacity 550 W at 15 $^{\circ}$ C (-10 25 $^{\circ}$ C)
- 2 F-308: cooling capacity 900 W at 15 °C (-10 25 °C)

Voltage

- 1 220 240 V
- 2 100 120 V
- ¹ PLASTIC+GLASS is sensitive to higher vapor temperatures
- $^{\rm 2}$ When ordering the glass assembly C, a recirculating chiller is not needed


Scope of delivery of Rotavapor® R-300 system configuration


All configurations are supplied ready to use.

Components	System configuration
Rotavapor® R-300 / R-300 Hand Lift**	1
Heating Bath B-301 / B-305**	1
Glass assembly	1
Vapor duct (with Combi-Clip)	1
Vacuum seal VS26	1
Evaporating flask (1 liter)	1
Receiving flask (1 liter)	1
Flyer "List of solvents"	1
Set of required cooling and vacuum tubings	1
Set of required power cables	1
Interface I-300 / I-300 Pro**	1
VacuBox	1
Communications cable. BUCHI COM 6 p	*
Vacuum Pump V-300	1
Recirculating Chiller F-305 / F-308**	*
Communication cables are included	*

^{*} Optional, according to order code

^{**} According to order code





Technical data

Rotavapor® R-300

	Rotavapor® R-300
Dimensions (W x D x H) with glass apparatus V	607 x 429 x 947 mm
Dimension (W x D x H) Shipping box	700 x 590 x 605 mm
Weight	13.5 kg (electric lift) 13.0 kg (manual lift)
Power supply	100 - 240 VAC ± 10%
Power consumption	100 W
Frequency	50/60 Hz
Ingress protection rating	IP21
Immersion angle	10 – 50°
Vertical adjustment range	220 mm
Limit position detection adjustment range	170 mm (electric lift) 100 mm (manual lift)
Rotation speed range	10 – 280 rpm
Max. flask capacity	3 kg
Approval	CE/CSA

Heating bath base, B-300

	Rotavapor® R-300
Dimensions (W x D x H)	183 x 286 x 82 mm
Weight	0.8 kg
Power supply	100 - 120 VAC ± 10% 220 - 240 VAC ± 10%
Frequency	50/60 Hz
Power consumption (with bath)	1500 W
Ingress protection rating	IP21
Approval	CE/CSA

Heating bath


	Heating Bath B-301	Heating Bath B-305
Dimensions (W x D x H)	225 x 212 x 192 mm	307 x 275 x 202 mm
Weight	2.0 kg	4.2 kg
Power supply	100 - 120 VAC ± 10% 220 - 240 VAC ± 10%	100 - 120 VAC ± 10% 220 - 240 VAC ± 10%
Frequency	50/60 Hz	50/60 Hz
Power consumption	1250 W	1500 W
Heating power	1100 W	1300 W

	Heating Bath B-301	Heating Bath B-305
Ingress protection rating	IP21	IP21
Temperature range	up to 95 °C	up to 220 °C
Evaporating flask maximum size	1000 mL	5000 mL
Maximum bath capacity	2.3 L	5.5 L
Adjustment accuracy	± 1 °C	± 1 °C
Regulation precision	At 60 °C: ±1 °C At 95 °C: ±1 °C	At 60 °C: ±1 °C At 95 °C: ±1 °C At 180 °C: ±3 °C At 220 °C: ±4 °C
Heating fluid	Water	Water/Oil
Approval	CE/CSA	CE/CSA

Heating bath performance

Heating-up time of Heating Bath B-301 and B-305.

	Water		Oil (Polyalkylenglyc	ole)
	50 °C	95 °C	180 °C	220 °C
B-301	4 min	13 min		
B-305	8 min	22 min	27 min	37 min



 $\Delta T_{70~^{\circ}\text{C (water)}}$: B-301 -> 13 min; B-305 -> 22 min

Special features

Rotavapor® R-300

'	
Firmware update	With the BUCHI Monitor App and BUCHI Bluetooth Dongle, firmware updates can be easily done by the user at no costs.
End stop positioner	Adjustable via button within a range of 170 mm (electronic) / 100 mm (manual)
	Prevents glass breakage and damages to heating bath pan
Multifunctional Combi- Clip	Tool for evaporating flask fixation and removal, as well as vapor duct release. Allows one handed flask removal
Immersion angle adjustable	Adjustable within 7-stages enables the use of different flask sizes, adjustable for individual application
Safety stop (electronic lift only)	Automatic operation stop by power failure Lifts flask automatically out of heating bath
2 lift speeds	Moves quickly until 2 cm above close to depth stop, it than moves slowly to depth stop Allows a smooth operation

Rotavapor® R-300

Earthquake safety	Lug on backside for bath and Rotavapor® Enables fixation on lab bench
Lift height indication	Position is digitally indicated on display of heating bath Allows reproducibility of parameter for future applications
Color code	Green cables indicate BUCHI communication cable, marked with "COM" Allows a simple understanding of set up and easy play & plug
Large top hole V-con- denser	Top hole with screw cap (SVL22) Allows easy cleaning of condenser
Ergonomic handle	Electronic and manual version available Allows to conveniently change the position of lift
Display	Rotation speed (actual value), lift position (actual value) and heating bath temperature (set and actual value) displayed on the Heating Bath B-300 Base. Allows user to monitor all parameters at a glance.

Rotavapor® related patents

Name	Date of application	Number of patent
Condensate trap	2016/06/29	Pending
Rotary evaporator for use in laboratory, has base for holding heat bath, where base includes electri- cal two-way plug-in connection unit which is de- tachably connected with plug-in connection unit of heat bath		DE 202005015235 U1
Height-adjustable support with holder and method of uncoupling the holder	1 2000/06/26	EP 1 484 112 EP 1 068 903 DE 1 068 903

Heating Bath B-301 / B-305

Universal use (B-305 only)	Useable for multiple applications up to 220 °C Allows the use of water and oil as a heating media
Overheat protection	Mechanical and electronic over-temperature control Cuts off power when actual exceeds set temperature and by rapid temperature increase
Carry handles	Two integrated laterally ergonomically designed carry handles Allows easy carrying for re-filling or emptying the heating bath, without the risk of burn
Bath guidance	Integrated in Heating Bath B-300 Base Allows an easy positioning of bath along horizontal axe, for individual application
Temperature limit fixation (B-305 only)	Prevents accidental change of set temperature (e.g. when water is used as heating media and maximal temperature is fixed at 95 °C)
Key-button locking functions	To lock the set temperature of the heating bath Prevents accidental change of set conditions during process
Infrared communication	Transition of heating bath data to interface Allows immaculate, fast communication

Kettle connection	Cordless bath
	Allows easy emptying and refilling of bath without unplugging
Heating Bath B-300 Base ex-	Functional with two different heating baths
tendable	Allows the use of B-301 or B-305, adequate to application
LCD-display	Digital, displays all operating parameters Indicates set and actual values of heating bath temperature, rotation speed and lift position

Optional accessories

Protection shield	Special shaped, robust shield for user protection enables handling on evaporating flask while being protected
Top cover for heating bath	To cover B-305 when not in use to safe energy and for less evaporation of heating medium, no drip off, convenient handle prevents burning skin

Compatibility and communication requirements with R-300

Interface / Vacuum Controller	Compatibility	Communication requirements
I-300 / I-300 Pro	compatible	VacuBox 2x Communication cable. BUCHI COM 6p: 11058705 (0.3 m), 11058707 (1.8 m), 11058708 (5 m) or 11064090 (15 m)
I-100	-	-
V-850 / V-855	-	-
V-800 / V-805	-	-

Explanation of terms

compatible	Rotavapor® and vacuum controller/interface communicate with each other.
-	Connection not possible

A VacuBox is included when ordering Rotavapor® R-300 with an Interface I-300 or Interface I-300 Procompatbility pack. Additional features for Rotavapor® R-300 are listed as follows.

Feature comparison


Features in combination with R-300	I-300 Pro	I-300	I-100	Old vacuum controllers (V-850 / V-855, V-800 / V-805)
Live data charting	•			
Digitally controlled lift position	•			
Data recording	•			
Guided operation (methods)	•			
Automatic ON/OFF rotation speed	•	•		
Automatic and manual up/down lift	•	•		
Display and setting rotation speed	•	•		

Features in combination with R-300	I-300 Pro	I-300	I-100	Old vacuum controllers (V-850 / V-855, V-800 / V-805)
Display and setting bath temperature	•	•		
Dynamic distillation (only when connected with Vacuum Pump and Recirculating Chiller)	•	•		
Remote monitoring with Monitor App	•	•		
ECO mode (shuts down Recirculating Chiller and Heating Bath)	•	•		
Automatic foam detection	•	•		
Management from a distance	•	•		

Accessories - Rotavapor® glass assemblies

The following items have to be purchased in addition to a system that is ordered without glass.

All glass assemblies include a 1 liter receiving flask (corresponding to order code - either without coating, with P+G or P+G-LT), required hoses and a ball joint clamp. Evaporating flask, vacuum seal VS26, vapor duct and condenser holder are not included. All condensers allow feeding of the evaporating flask via stopcock.




429	607 630 750 608	615 593	808 607
Name	Type / Application	Cooling Surface	Order no.
		area (max.)	
Vertical: V	Vertical condenser	1500 cm ²	11062432
7777	 For standard applications, the most common condenser Connection for vapor temperature sensor Automatic distillation with Interface I-300 / I-300 Pro possible Use of foam sensor possible 		11062433 P+G
Diagonal: A	Diagonal condenser	1500 cm ²	048168
F 2	For standard applications where height is limited		048169 P+G
Cold trap: C	Dry ice condenser	500 cm ²	040640
	 For distillation of solvents with low boiling points Maximum condensation of vapors No cooling water needed, but i.e. dry ice or ice Use of foam sensor possible 		040642 P+G
Reflux: S	Vertical reflux condenser with shut-off valve	1500 cm ²	048290
	 For combining reflux reaction and distillation Connection for vapor temperature sensor Automatic distillation with Interface I-300 / I-300 Pro possible Use of foam sensor not possible 		048291 P+G- LT

Name	Type / Application	Cooling Surface	Order no.
		area (max.)	
Cold trap reflux: CR	 Dry ice condenser with shut-off valve For combining reflux reaction and distillation For distillation of solvents with low boiling points Maximum condensation of vapors Connection for vapor temperature sensor No cooling water needed, but i.e. dry ice or ice Use of foam sensor not possible 	500 cm ²	048292 048293 P+G- LT
Expansion: E	Descending condenser with expansion vessel	1500 cm ²	11061112
965	 Ideal for distillations exhibiting foaming and "bumping" tendencies/products Connection for vapor temperature sensor Automatic distillation with Interface I-300 / I-300 Pro possible 		11061113 P+G
Double	Vertical reflux condenser with double jacket for	1300 cm ²	048176
jacket: BY	 Intensive condensation and with shut-off valve Vertical intensive condenser with double jacket and distributor for particularly efficient condensation For combining reflux reaction and distillation For particularly efficient and intensive condensation Connection for vapor temperature sensor 		048297 P+G
	 With an additional SJ 29/32 joint on top for flexi- ble expansion 		
	Use of foam sensor not possible		
High performance: HP	 High performance condenser for increased distillation rate Double cooling capacity as standard V condenser Enables distillation at lower pressure values and higher temperature differences Reduces emissions of solvents Reduces the wear of vacuum pump parts Compatible with: AutoDest sensor Foam sensor Vapor temperature sensor Condensate trap 	3000 cm ²	11066561 11066562 P+G

For more information concerning the glass assemblies please check the "Laboratory Evaporation Glassware" brochure.

Features of the glass assemblies



- 1 Top inletFor easy cleaning of the condenser
- 2 Inner glass tube For reduced chimney effect upon vacuum applications
- 3 Stopcock (available for R-300, R-215, R-210, R-100, R-3; see FAQ)
 - Used for aeration and feeding (by using the feeding hose)
- 4 Inner glass plate
 Prevention of solvent back flow
- 5 Connector for sensor
 For measurement of vapor temperature during distillation
- 6 Condensate trap

Characteristics	V	С	Α	CR	S	Е	BY	HP
Top inlet	• (SVL22; extra large)		•		•	•	•	• (SVL22; extra large)
Inner glass tube	•				•		·	•
Flask feeding via stopcock possible	•	•	•	•	•	•	•	•
Inner glass plate	•							•
Connector for va- por temperature sensor	•				•	•	•	•
Condenser holder	01	01	01	01	01	•	01	• (048180)
Reflux possibility (part included)				•	•		•	
Condensate trap	•							•

o is optional

is included

¹ Condenser holder cpl. (048180) is recommended

Accessories

Heating baths

	Order no.	Image
Heating Bath B-301, 110V, not incl. Base B-300 For up to 1 Liter evaporating flasks. For temperatures up to 95°C. Meant to be used with water as heating media.	11B301002	
Heating Bath B-301, 230V, not incl. Base B-300 For up to 1 Liter evaporating flasks. For temperatures up to 95 °C. Meant to be used with water as heating media.	11B301001	
Heating Bath B-301, 110V, incl. Base B-300 For up to 1 Liter evaporating flasks. For temperatures up to 95 °C. Meant to be used with water as heating media.	11B301102	
Heating Bath B-301, 230V, incl. Base B-300 For up to 1 Liter evaporating flasks. For temperatures up to 95 °C. Meant to be used with water as heating media.	11B301101	
Heating Bath B-305, 110V, not incl. Base B-300, 110V For up to 5 Liter evaporating flasks. For temperatures up to 220 °C. Meant to be used with water or oil as heating media.	11B305002	
Heating Bath B-305, 230V, not incl. Base B-300 For up to 5 Liter evaporating flasks. For temperatures up to 220 °C. Meant to be used with water or oil as heating media.	11B305001	
Heating Bath B-305, 110V, incl. Base B-300 For up to 5 Liter evaporating flasks. For temperatures up to 220 °C. Meant to be used with water or oil as heating media.	11B305102	
Heating Bath B-305, 230V, incl. Base B-300 For up to 5 Liter evaporating flasks. For temperatures up to 220 °C. Meant to be used with water or oil as heating media.	11B305101	

	Order no.	Image
Heating Bath B-305, 110V, for R-21x, incl. Base	11B305112	
For up to 5 Liter evaporating flasks. For temperatures up to 220 °C. Meant to be used with water or oil as heating media. Included: Adapter for Rotavapor® R-21x.		
Heating Bath B-305, 230V, for R-21x, incl. Base	11B305111	
For up to 5 Liter evaporating flasks. For temperatures up to 220 °C. Meant to be used with water or oil as heating media. Included: Adapter for Rotavapor® R-21x.		

Heating bath accessories

riodding bath accessories		
	Order no.	Image
Adapter B-300. For Heating Bath B-305	11061317	_
Required when Heating Bath B-305 is used with a Rotavapor® R-210, R-215 or R II. Not compatible with Heating Bath B-301.		
Dewar accessory, set. Incl. Dewar container and adapter	11066662	
For sample preparation in freeze drying. To be used with dry ice and ethanol / isopropanol / acetone or separately with liquid nitrogen. Compatible with Heating Bath B-300 (not plugged in during operation) and Rotavapor® R-300.		
Dewar container, not incl. adapter	11066645	
For sample preparation in freeze drying. To be used with dry ice and ethanol / isopropanol / acetone or separately with liquid nitrogen. To be used only with adapter, Heating Bath Base B-300 (not plugged-in during operation) and Rotavapor R-300.		
Lid. For Heating Bath B-305	11059500	<u> </u>
To save energy and to minimize loss of water when heating bath is not in use.		
Protection shield, Operator protection for heating bath	11061402	
For Heating Bath B-301 and B-305. Material: 4 mm thick PMMA.		
Splash shield, Operator protection against splashes	11064768	
For Heating Bath B-301 and B-305. Material: 6 mm thick polycarbonate.		

Protective cover. For Heating Bath B-300 base

11064111

Used for protection of the heating bath display.



Further accessories heating baths

	Order no.
Heating bath balls. 450 pcs., PP, Ø10 mm	036405
To reduce energy consumption of heating bath and for less evaporation of the heating	
medium. For temperatures up to 100 °C.	

Tubing

	Order no.
Tubing. FEP, Ø6/8 mm, transparent, per m	027900
Use: Vacuum, cooling media.	
Tubing. Natural rubber, Ø6/16 mm, red, per m	017622
Use: Vacuum.	
Tubing. Nylflex, PVC-P, Ø8/14 mm, transparent, per m	004113
Tubing. PTFE, Ø4.7/5.5 mm, transparent, 330 mm	000646
Use: To introduce solvent into evaporating flask during distillation.	
Tubing. PTFE, Ø4.7/5.5 mm, transparent, 460 mm	000643
Use: To introduce solvent into evaporating flask during distillation.	
Tubing. PTFE, Ø8/10 mm, white, per m	027277
Use: Vacuum, feeding (industrial Rotavapor®).	
Tubing. Silicone, Ø6/9 mm, transparent, per m	004133
Use: Cooling media.	
Tubing. Synthetic rubber, Ø6/13 mm, black, per m	11063244
Use: Vacuum.	

Insulation tubing

	Order no.	Image
Hose insulation, Kaiflex, 11/23 mm, 1 m, black	028696	

Further accessories

	Order no.	Image
StatusLight cpl., incl. communication cable	11068959	
Indicates the status of the instrument (instrument is ready to use, has an error or is in operation).		

	Order no.	Image
BUCHI Bluetooth® Dongle, connects instrument to smart- phone via Bluetooth®	11067770	
Needed for firmware updates and eSupport.		
AutoDest sensor. Incl. cap nut, seal GL14	11059225	
For automatic distillation. Measures temperature of cooling media and the vapor temperature. Vacuum is adjusted according to cooling capacity of condenser. Meant to be used with the Interface I-300 / I-300 Pro and glass assembly V, HP or S.		
Condensate trap.	11062955	
Condenser V and HP, grey, TPE, hose connector Ø8 mm. Collects and drains condensate which may accumulate on the condenser.		
Condensate trap.	11061985	
Condenser V and HP, green, TPE, hose connector Ø8 mm. Collects and drains condensate which may accumulate on the condenser.		
Cooling water valve. 24VAC	031356	~ m
Valve opens cooling water feed during distillation. Meant to be used with a vacuum controller/interface.		
Extraction thimbles, Set. 25 pcs, 43 x 118 mm, cellulose For Soxhlet extraction unit.	018106	
Flask holder. EPDM, slip free	048618	
Holder for round-bottom flasks (50-5000 mL).	0.00.0	The state of the s
Flask holders, set. 5 pcs., EPDM, slip free	11059916	
Holder for round-bottom flasks (50-5000 mL).		
Foam sensor. Incl. holder	11061167	
Prevents sample from foaming into the condenser by aeration of system temporarily. Meant to be used with the Interface I-300 / I-300 Pro.		

	Order no.	Image
IQ/OQ. R-300 en	11063608	(1) The State of State (1)
Official BUCHI documents.		
Lab grease. Glisseal 40 (60 g)	048197	
To grease joints in order to increase tightness of system.		
Level sensor. Incl. O-ring, cable, rubber band	11060954	
Prevents an overflow in the receiving flask of condenser/secondary condenser. Placed on receiving flask. Meant to be used with the Interface I-300 / I-300 Pro.		
Repeating OQ. R-300 en	11063609	Relationary A. Salah
Official BUCHI documents.		
Power plug fixation. For R-300	11063710	
Used to fixate the power plug into Rotavapor®.		
Set earthquake fastener. Lug for backside	11062386	
To fasten device on lab bench.		eoe
Thread adapter. For vacuum connection, SVL22/GL14	040871	
Use: To switch vacuum connection from a Rotavapor® to Syncore® vacuum distributor (037686).		
Vacuum connection. Incl. drain valve, hose barb Ø 8 mm, BJ35/20	001006	
Connection piece for aeration of the system, placed between condenser and receiving flask.		
Vacuum valve. Magnetic valve, 24V/4W, Mini-DIN, 1.5 m	11060706	
Flow valve without flask, meant to be used with a centralized vacuum source or an unregulated vacuum pump. Meant to be used with the Interface I-300 /I -300 Pro.		

	Order no.	lmage
Vapor temperature sensor. Incl. cap nut, seal GL14	11060707	
Measures the vapor temperature inside the system. Meant to be used with the Interface I-300 / I-300 Pro.		
Water jet pump. Plastic	002913	
Used when tap water is used to generate vacuum.		
Water regulation nozzle. Flow regulator, incl. hose clamp, sieve	011606	
Used when tap water is used to generate vacuum. Reduces water consumption.		
Woulff bottle. 125 mL, P+G, incl. holder	11059905	P —
For trapping particles and droplets and for pressure equalization.		
Woulff bottle. 3-neck, 800 mL, P+G.	025519	
For trapping particles and droplets and for pressure equalization.		
Manometer with needle valve. Incl. vacuum gauge, not incl. holder	047391	
Accessories for glass assemblies		
	Order no.	Image
Condenser holder.For glass assembly V/HP/C/S/CR/BY	048180	
Content: Holding rod, rubber band, cross sleve.		
Condenser holder. For Glass assembly E	048181	
Content: Holding rod, bride, cross sleeve.		

	Order no.	Image
Extraction unit Soxhlet, 200 mL, incl. extraction part, reduction part	011744	
For Soxhlet extraction applications with a Rotavapor® R-300. Meant to be used with glass assembly S. Condenser holder is included.		
Extraction unit Soxhlet, 500mL, incl. extraction part, reduction part	011745	
For Soxhlet extraction applications with a Rotavapor® R-300. Meant to be used with glass assembly S. Condenser holder is included.		

Stopcocks

	Order no.	Image
For condenser C/CR, glass, SJ18.8/38	040628	
For aeration of the system. For cold trap outer part.		
Standard, glass, SJ18.8/38	040627	
For aeration of the system.		
Stopcock, Analytic PTFE/25% glass fiber, SJ18.8/38	11069607	
For feeding of solvents and aerating the system.Less cross-contamination compared to standardstopcock. For applications where grease should be avoided.		
Content: PTFE stopcock (no tubing included).		U
PTFE, incl. 3-way valve	11058814	
For feeding of solvents and aeration of the system. For applications when grease should be avoided. Used instead of standard-stopcock (040627). Content: Inlet tubing 300mm, backfeed tubing 600mm, cap nut GL10.		

Wearing parts

Vacuum seal

	Order no.
Vacuum seal VS26, PTFE base, NBR O-ring, FDA-complaint	11069167

Seals

	Order no.
Seal	038225
For cap nut GL14, EPDM seal	
Seals, set. 10 pcs, for hose barbs GL14, EPDM, black	040029

	Order no.
Seals, set. 10 pcs, for hose barbs GL14, FPM, green	040040
Seals, set. 10 pcs, for hose barbs GL14, silicone, red	040023

Hose barbs

Hose barbs		
	Order no.	Image
Hose barb. Bent, GL14, incl. silicone seal	018916	
Hose barbs, set. 2 pcs, bent (1), straight (1), GL14, silicone seal Content: Hose barbs, cap nuts, seals	041939	
Hose barbs. set. 3 pcs, bent, GL14, silicone seal Content: Hose barbs, seals.	041987	
Hose barbs, set. 4 pcs, bent GL14, silicone seal Content: Hose barbs, cap nuts, seals	037287	O COLOR OF THE PROPERTY OF THE
Hose barbs. set. 4 pcs, bent, GL14, EPDM seal Content: Hose barbs, cap nuts, seals.	043129	TO THE PARTY OF TH
Hose barbs, set 4 pcs, bent, GL14, FPM seal Content: Hose barbs, cap nuts, seals.	040295	
Hose barbs. set. 4 pcs, straight, GL14, EPDM seal Content: Hose barbs, cap nuts, seals.	043128	
Hose barbs. set. 4 pcs, straight, GL14, FPM seal Content: Hose barbs, cap nuts, seals.	040296	

	Order no.	lmage
Hose barbs, set. 4 pcs, straight, GL14, silicone seal	037642	iriago
Content: Hose barbs, cap nuts, seals		Company of the compan
Hose barbs. set. 6 pcs, bent (4), straight (2), GL14, silicone seal Content: Hose barbs, cap nuts, seals.	038000	O GO
Miscellaneous wearing parts	Order no.	lmage
Cap nut. Screw cap with hole GL10	023875	
Cap nuts, set. 10 pcs, screw cap with hole GL14, seal EPDM Content: Cap nuts, seals	041999	
Cap nuts. set. 10 pcs, screw cap with hole, GL14	041956	
Hose, set. Incl. GL14 set, FEP tube 11061756	11065373	
Used to connect VacuBox and Valve unit/Woulff bottle when both are fixed on the Rotavapor® R-300.	11000070	
O-ring. FKM, 40.9/2.6 mm, for reflux insert	048078	
O-ring. FPM, for cap nut GL10, Ø3.00/2.70 mm	023900	

	Order no.	Image
Screw caps. set. 5 pcs, GL14	040624	
Screwing set, Hose barbs and seals GL14	11061921	
Content: Hose barbs bent 4 pcs, hose barbs straight 1 pc, cap nuts 4 pcs, screw caps 2 pcs, seals EPDM 3 pcs.		O O O O O O O O O O O O O O O O O O O

Rotavapor® glass accessories - Vapor ducts









Vapor ducts compatible with VS26 (incl. Combi-Clip)*

SJ	29/32	24/40	29/42	24/29
For glass assembly V, C, S, E, CR, BY, HP	11062186	11062187	11062464	11062909
For glass assembly A	11062267	11062268	11062269	
For glass assembly V, C, HP (analytical)	11062465	11062466	11062467	
Vapor duct with frit, Ø 26 mm, SJ29/32, incl. Combi-Clip	11057297			
For powder drying. To prevent powder from getting into the condenser. For glass assembly V, C, S, E, BY, HP and CR.				
For high temperatures 29/32, short cpl.	11061837			
Content: Combi clip, vapor duct HT, ring NS 34/32 on 30/32				

^{*}Single Combi-Clip: 11059770

Additional glassware

Additional glassware like evaporating flasks, drying flasks, receiving flasks, beaker flasks, Reitmeyer adapters (bump traps) and distillation spiders can be found in the Laboratory Evaporation Glassware brochure.