855 Robotic Titrosampler



Manual 8.855.8001EN / 2020-02-29





Metrohm AG CH-9100 Herisau Switzerland Phone +41 71 353 85 85 Fax +41 71 353 89 01 info@metrohm.com www.metrohm.com

855 Robotic Titrosampler

Manual

8.855.8001EN / 2020-02-29

Technical Communication Metrohm AG CH-9100 Herisau techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Table of contents

1	Introd	uctio	n	1
		1.1	The 855 Robotic Titrosampler in the Titrando system	. 1
		1.2 1.2.1 1.2.2 1.2.3	Instrument description Model versions Instrument components Intended use	3 4
		1.3 1.3.1	About the documentation	
		1.4 1.4.1 1.4.2 1.4.3 1.4.4 1.4.5 1.4.6	Safety instructions General notes on safety Electrical safety Tubing and capillary connections Personnel safety Flammable solvents and chemicals Recycling and disposal	6 6 7 7 9
2	Overv	iew o	f the instrument	10
		2.1	Front and rear	10
		2.2	Back panel	12
		2.3	Sensor connectors	12
		2.4	Sample racks	13
3	Install	ation		15
		3.1 3.1.1 3.1.2 3.1.3	Setting up the instrument Packaging Checks Location	15 15
		3.2 3.2.1	Preparing the Sample Processor Connecting the instrument to the power grid	
		3.3	Connecting a computer	16
		3.4	Connecting the Swing Head	18
		3.5	Configuring the robotic arm	19
		3.6	Mounting the robotic arm	21
		3.7	Robotic arms with beaker sensor	23
		3.8	Installing the rinsing and aspiration equipment	24
		3.9	Placing tubing and cables in the guide chain	27
		3.10	Equipping the titration head	28

		3.11	Connecting the tower stirrer	31
		3.12	Connecting an external pump	32
		3.13 3.13.1 3.13.2 3.13.3	Connecting MSB devices Connecting a dosing device Connecting a stirrer or titration stand Connecting a Remote Box	34 35
		3.14 3.14.1	Connecting USB devices Connecting a barcode reader	
		3.15	Mounting the base plate	39
		3.16	Mounting the drip pan	40
		3.17	Attaching a sample rack	41
		3.18	Mounting the safety shield	42
4	Opera	ation a	nd maintenance	44
		4.1	General notes	44
		4.2	Care	44
5	Troub	lesho	oting	45
		5.1	Sample Processor	45
		5.2	Robotic arm	45
		5.3	Pump	46
6	Appe	ndix		47
		6.1	Beaker sensor	47
		6.2	Rinsing nozzles	47
		6.3 6.3.1	Remote interface Pin assignment of the remote interfaces	
		6.4	Stirring rate	50
		6.5 6.5.1 6.5.2 6.5.3	Robotic arms Robotic arms for titration Robotic arms for sample preparation Robotic arms for special applications	. 51 53
7	Techn	ical sp	pecifications	55
		7.1 7.1.1 7.1.2 7.1.3	Measuring interface Potentiometry Polarizer Temperature	. 55 55
		7.2	Lift and turntable	56
		7.3	Membrane pump(s) with valve	57

	7.4	786 Swing Head	57
	7.5	Interfaces and connectors	57
	7.6	Power connection	58
	7.7	Ambient temperature	58
	7.8	Reference conditions	58
	7.9	Dimensions	58
8	Accessories	5	59
	Index		60

Table of figures

Figure 1	The Titrando system	1
Figure 2	Front 855 Robotic Titrosampler	10
Figure 3	Rear 855 Robotic Titrosampler	11
Figure 4	Connector strip	12
Figure 5	Sensor connectors	13
Figure 6	Sample rack 6.2041.840	13
Figure 7	Connecting the computer	17
Figure 8	Connecting the Swing Head	18
Figure 9	Robotic arms - standard model versions:	19
Figure 10	Configuration data of the robotic arms	20
Figure 11	Limitation screw on the robotic arm	21
Figure 12	Mounting the robotic arm	
Figure 13	Connecting a beaker sensor (for example 6.1462.150)	23
Figure 14	Mounting the rinsing and aspiration tubings	24
Figure 15	Installing the rinsing tubings and the aspiration tip	29
Figure 16	Installing stirrer and electrode	
Figure 17	Rod stirrer 802 Stirrer	31
Figure 18	Magnetic stirrer 741 Stirrer	31
Figure 19	Connecting the tower stirrer	31
Figure 20	Connecting a pump	32
Figure 21	MSB connections	33
Figure 22	Connecting a dosing device	35
Figure 23	Connecting an MSB stirrer	
Figure 24	Connecting the rod stirrer to the titration stand	36
Figure 25	Connecting the Remote Box	37
Figure 26	USB connectors	
Figure 27	Mounting the base plate	39
Figure 28	Installing the drip pan	
Figure 29	Attaching rack	41
Figure 30	Mounting the safety shield	42
Figure 31	Beaker sensor on the tower	
Figure 32	Spray nozzles - Mode of operation	
Figure 33	Connectors of the Remote Box	
Figure 34	Pin assignment of remote socket and remote plug	
Figure 35	Rotational speed depending on the stirring rate	50

1 Introduction

1.1 The 855 Robotic Titrosampler in the Titrando system

The 855 Robotic Titrosampler is a component of the modular Titrando system. Operation is carried out either by a Touch Control with a touch-sensitive screen ("stand-alone titrator") or by a computer with a corresponding software.

A Titrando system can contain numerous, various kinds of instruments. The following figure provides an overview of the peripheral devices you can connect to the 855 Robotic Titrosampler.



Up to three control instruments (Titrando, Dosing Interface, USB Sample Processor, etc.) can be controlled via USB connection during operation with the 900 Touch Control. If the tiamo software is used, then there is no limit to the number of control instruments that can be used to extend the system.

Updating the device software is described in the manual for the Touch Control or in the online help for the tiamo software.

1.2 Instrument description

The 855 Robotic Titrosampler is a combined analyzer, unifying the functionality of a titrator and a sample changer. It can be seamlessly incorporated as a system component and control instrument into a Metrohm Titrando system.

The 855 Robotic Titrosampler is a versatile instrument. It has been designed exclusively for usage in factories, as long as ambient conditions *(see chapter 7.7, page 58)* are fulfilled, and laboratories and thereby covers a wide range of applications.

Thanks to the integration of high-performance USB interfaces, it can be incorporated seamlessly into a Metrohm Titrando system. The versatile communications options of the Titrando system (Remote Box, LIMS link, etc.) can thus all be used. Thanks to these abilities, a 855 Robotic Titrosampler is predetermined for all kind of automation tasks in a modern laboratory, especially for highly integrated laboratory data systems.

The user interface of the Touch Control or the *tiamo* software guarantees comfortable operation and programming of the 855 Robotic Titrosampler. The comprehensive range of commands and the various configuration possibilities can comfortably and efficiently be used this way. The integration into the Titrando system also guarantees a 100% conformity of the complete automation system with the FDA regulations (Food and Drug Administration), especially with the regulation 21 CFR part 11, electronic records and signatures.

There are exchangeable standard sample racks available for many vessel dimensions. Freely selectable special beaker positions can be defined for e.g. rinsing or conditioning beakers on every rack.

The equipment with a 786 Swing Head allows to process a large number of samples in one single sample series. The robotic arms for the 786 Swing Head make it possible to move to any given point on a sample rack. This way the number (a maximum of 999 rack positions) and sequencing of the samples is almost unlimited.

Customer-specific special racks for individual requirements can be fabricated upon request.

1.2.1 Model versions

The 855 Robotic Titrosampler is available in the following versions with different components.

2.855.0010

1-tower version

- 1 membrane pumpe and 1 valve
- 1 connector for an external pump
- 1 stirrer connector (tower stirrer)
- 1 786 Swing Head
- 1 measuring interface
- 3 MSB connectors for dosing devices, Remote Box or stirrers
- 2 USB connectors
- 1 connector "Controller" for Touch Control or PC

2.855.0020

1-tower version

- 2 membrane pumps and 2 valves
- 1 stirrer connector (tower stirrer)
 - 1 786 Swing Head
 - 1 measuring interface
- 3 MSB connectors for dosing devices, Remote Box or stirrers
- 2 USB connectors
- 1 connector "Controller" for Touch Control or PC

2.855.0030

1-tower version

• 1 membrane pumpe and 1 valve

- 1 connector for an external pump
- 1 stirrer connector (tower stirrer)
- 1 Swing Head connector
- 1 measuring interface
- 3 MSB connectors for dosing devices, Remote Box or stirrers
- 2 USB connectors
- 1 connector "Controller" for Touch Control or PC







2.855.0040 1-tower version 2 membrane pumps and 2 valves 1 stirrer connector (tower stirrer) 1 Swing Head connector

- 1 measuring interface
- 3 MSB connectors for dosing devices, Remote Box or stirrers
- 2 USB connectors
- 1 connector "Controller" for Touch Control or PC

1.2.2 Instrument components

The 855 Robotic Titrosampler has the following components:

Turntable

For sample racks with a diameter of up to 48 cm.

One tower with lift

The 2.855.0010 and 2.855.0020 versions with mounted 786 Swing Head. The 2.855.0030 and 2.855.0040 versions can be upgraded with a Swing Head. Any robotic arm can be mounted.

- One or two membrane pumps per tower
 Instead of an integrated pump, an external pump connector is available depending on the model version.
- One stirrer connector on the tower

For connecting a rod stirrer (802 Stirrer) or a magnetic stirrer (741 Stirrer).

Sensor connectors

One connector each for:

- potentiometric electrodes (pH, ISE, metal)
- reference electrodes
- polarizable electrodes
- temperature sensors (Pt1000 or NTC)
- Two USB connectors

For connecting a printer, barcode reader or other control instruments (Titrando, Dosing Interface etc.).

• Three MSB connectors (Metrohm Serial Bus)

For connecting dosing devices (Dosimat with exchange unit or Dosino with dosing unit), stirrers or Remote Boxes.

1.2.3 Intended use

The 855 Robotic Titrosampler is designed for usage as an automation system in analytical laboratories or in factories that fulfill the specified ambient conditions (*see chapter 7.7, page 58*). It is **not** suitable for usage in biochemical, biological or medical environments in its basic equipment version.



This instrument is suitable for processing chemicals and flammable samples. Usage of the 855 Robotic Titrosampler therefore requires the user to have basic knowledge and experience in handling toxic and caustic substances. Knowledge with respect to the application of the fire prevention measures prescribed for laboratories is also mandatory.

1.3 About the documentation

CAUTION

Please read through this documentation carefully before putting the instrument into operation. The documentation contains information and warnings which the user must follow in order to ensure safe operation of the instrument.

1.3.1 Symbols and conventions

The following symbols and formatting may appear in this documentation:

(5- 12)	Cross-reference to figure legend	
	The first number refers to the figure number, the sec- ond to the instrument part in the figure.	
1	Instruction step	
	Carry out these steps in the sequence shown.	
Method	Dialog text, parameter in the software	
File ► New	Menu or menu item	
[Next]	Button or key	
$\widehat{\mathbf{A}}$	WARNING	
	This symbol draws attention to a possible life-threat- ening hazard or risk of injury.	
	WARNING	
	This symbol draws attention to a possible hazard due to electrical current.	
	WARNING	
	This symbol draws attention to a possible hazard due to heat or hot instrument parts.	

	WARNING
	This symbol draws attention to a possible biological hazard.
	CAUTION
	This symbol draws attention to possible damage to instruments or instrument parts.
-	NOTE
	This symbol highlights additional information and tips.

1.4 Safety instructions

1.4.1 General notes on safety



WARNING

Operate this instrument only according to the information contained in this documentation.

This instrument left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.

1.4.2 Electrical safety

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.



WARNING

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.



WARNING

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

Supply voltage



WARNING

An incorrect supply voltage can damage the instrument.

Only operate this instrument with a supply voltage specified for it (see rear panel of the instrument).

Protection against electrostatic charges



WARNING

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Do not fail to pull the power cord out of the power socket before you set up or disconnect electrical plug connections at the rear of the instrument.

1.4.3 Tubing and capillary connections



CAUTION

Leaks in tubing and capillary connections are a safety risk. Tighten all connections well by hand. Avoid applying excessive force to tubing connections. Damaged tubing ends lead to leakage. Appropriate tools can be used to loosen connections.

Check the connections regularly for leakage. If the instrument is used mainly in unattended operation, then weekly inspections are mandatory.

1.4.4 Personnel safety



WARNING

Wear protective goggles and working clothes suitable for laboratory work while operating the 855 Robotic Titrosampler. It is also advisable to wear gloves when caustic liquids are used or in situations where glass vessels could break.



WARNING

Always install the safety shield supplied with the equipment before using the instrument for the first time. Pre-installed safety shields are not allowed to be removed.

The 855 Robotic Titrosampler may not be operated without a safety shield!



WARNING

Personnel are not permitted to reach into the working area of the instrument while operations are running!

A considerable risk of injury exists for the user.



WARNING

In the event of a possible blockage of a drive, the power plug must be pulled out of the socket immediately. Do not attempt to free jammed sample vessels or other parts while the device is switched on. Blockages can only be cleared when the instrument is in a voltage-free status; this action generally involves a **considerable risk of injury**.



WARNING

The 855 Robotic Titrosampler is **not** suitable for utilization in biochemical, biological or medical environments in its basic equipment version.

Appropriate protective measures must be implemented in the event that potentially infectious samples or reagents are being processed.

1.4.5 Flammable solvents and chemicals



WARNING

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set up the instrument in a well-ventilated location (e.g. fume cupboard).
- Keep all sources of flame far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

1.4.6 Recycling and disposal



This product is covered by European Directive 2012/19/EU, WEEE – Waste Electrical and Electronic Equipment.

The correct disposal of your old instrument will help to prevent negative effects on the environment and public health.

More details about the disposal of your old instrument can be obtained from your local authorities, from waste disposal companies or from your local dealer.

2 Overview of the instrument

2.1 Front and rear



Figure 2 Front 855 Robotic Titrosampler

- 1 Safety shield (6.2751.100) Other models, *see chap. Accessories*.
- **3 Turntable** With guide bolts.
- 5 Guide chain For cables and tubings.
- 7 786 Swing Head (2.786.0020) Drive for the robotic arm.
- 9 Beaker sensor

- 2 Sample rack (6.2041.800) Other models, *see chap. Accessories*.
- 4 Assembly rail For magnetic stirrer (741 Stirrer).
- 6 Lift With titration head holder.
- 8 Robotic arm (6.1462.050) With aspiration and dosing tips.



Figure 3 Rear 855 Robotic Titrosampler

- Distributor For rinsing equipment.
 Pump connector
 - Pump 2. For external pumps (e.g. 772 Pump Unit or 823 Pump Unit).

7

- 5 Connector "Swing Head" Mini DIN socket (9-pin).
- **7** Back panel with connectors
- **9 Warning symbols** See chapter *Safety instructions*.

- 2 Membrane pump Pump 1.
- 4 Pump valve
- 6 Stirrer connector DIN socket. For rod stirrer (802 Stirrer) or magnetic stirrer (741 Stirrer).
- 8 Measuring interface With various measuring inputs.

2.2 Back panel



1USB connector (USB 1 and USB 2)23MSB connector (MSB 1 to MSB 3)
For stirrers, dosing devices, Remote Box.45Type plate
Contains specifications concerning supply
voltage and serial number.5

2 Connector (Controller)

For the connection to the PC or Touch Control.

4 Power socket

2.3 Sensor connectors

The connectors for electrodes and temperature sensors are on the rear of the tower.



Figure 5 Sensor connectors

- 1 Electrode connector (Ref.) For connecting reference electrodes, e.g. Ag/AgCl reference electrode. Socket B, 4 mm.
- **3** Temperature sensor connector (Temp.) For connecting temperature sensors (Pt1000 or NTC). Two B sockets, 2 mm.
- 2 Electrode connector (Ind.) For connecting pH, metal or ion-selective electrodes with integrated or separate reference electrode. Socket F.
- 4 Electrode connector (Pol.) For connecting polarizable electrodes, e.g. double Pt electrodes. Socket F.

2.4 Sample racks

A sample rack is a turntable that acts as a receptacle for sample vessels. Various types of sample racks are available for different numbers and types of sample vessels.

The 855 Robotic Titrosampler requires sample racks with up to a maximum of **48 cm diameter**.



Figure 6 Sample rack 6.2041.840

Other user-defined sample racks can be supplied upon request and the required rack data can be loaded and configured in the control software. Any arrangement of rack positions is possible.

Magnet codes

Every single sample rack can be unambiguously identified by means of a magnet code. The Sample Processor can thus recognize automatically which rack is attached.

When replacing a sample rack, this should first be returned to starting position using the **[Reset rack]** function (see "Manual Operation" in the control software).

This will enable an unambiguous recognition of the rack and thus the correct positioning of the beaker. A positioning table is assigned to each rack type in which each rack position is defined.

3 Installation

3.1 Setting up the instrument

3.1.1 Packaging

The instrument is supplied in protective packaging together with the separately packed accessories. Keep this packaging, as only this ensures safe transportation of the instrument.

3.1.2 Checks

Immediately after receipt, check whether the shipment has arrived complete and without damage by comparing it with the delivery note.

3.1.3 Location

The instrument has been developed for operation indoors and may not be used in explosive environments.

Place the instrument in a location of the laboratory which is suitable for operation and free of vibrations and which provides protection against corrosive atmosphere and contamination by chemicals.

The instrument should be protected against excessive temperature fluctuations and direct sunlight.

3.2 Preparing the Sample Processor

3.2.1 Connecting the instrument to the power grid



WARNING

Electric shock from electrical potential

Risk of injury by touching live components or through moisture on live parts.

- Never open the housing of the instrument while the power cord is still connected.
- Protect live parts (e.g. power supply unit, power cord, connection sockets) against moisture.
- Unplug the power plug immediately if you suspect that moisture has gotten inside the instrument.
- Only personnel who have been issued Metrohm qualifications may perform service and repair work on electrical and electronic parts.

Connecting the power cord

Accessories

- Power cord with the following specifications:
- Length: max. 2 m
- Number of cores: 3, with protective conductor
- Instrument plug: IEC 60320 type C13
- Conductor cross-section 3x min. 0.75 mm² / 18 AWG
- Power plug:
 - according to customer requirement (6.2122.XX0)
 - min. 10 A



NOTICE

Do not use a not permitted power cord!

1 Plugging in the power cord

- Plug the power cord into the instrument's power socket.
- Connect the power cord to the power grid.

3.3 Connecting a computer

The 855 Robotic Titrosampler requires a USB connection to a computer in order to be able to be controlled by a computer software. Using a 6.2151.000 controller cable, the instrument can be connected directly, either to a USB socket on a computer, to a connected USB hub or to a different Metrohm control instrument.

You need administrator rights for the installation of driver software and control software on your computer.

Cable connection and driver installation

A driver installation is required in order to ensure that the 855 Robotic Titrosampler is recognized by the computer software. To accomplish this, you must comply with the procedures specified. The following steps are necessary:

1 Install the software

- Insert the computer software installation CD and carry out the installation program directions.
- Exit the program if you have started it after the installation.

2 Establishing the cable connections

- Connect all peripheral devices to the instrument, *see chapter* 3.13, page 33.
- Connect the instrument to the power grid if you have not already done this (see chapter 3.2.1, page 15).
- Connect the instrument to a USB connector (type A) of your computer (see manual of your computer). The 6.2151.000 cable is used for this purpose.



Figure 7 Connecting the computer

The instrument is recognized. Depending on the version of the Windows operating system used, the driver installation proceeds differently afterwards. Either the necessary driver software is installed automatically or an installation wizard is started.

3 Follow the instructions of the installation wizard.

If problems should occur during installation, contact your company's IT support team.



NOTICE

The plug on the instrument end of the 6.2151.000 controller cable is protected against accidental disconnection by means of a pull-out protection feature. If you wish to pull out the plug, you first need to pull back the outer plug sleeve marked with arrows.



Registering and configuring the instrument in the computer software

The instrument must be registered in the configuration of your computer software. Once that has been done, you can then configure it according to your requirements. Proceed as follows:

1 Setting up the instrument

- Start the computer software. The instrument is automatically recognized. The configuration dialog for the instrument is displayed.
- Make configuration settings for the instrument and its connectors.

More detailed information concerning the configuration of the instrument can be found in the documentation for the respective computer software.

3.4 **Connecting the Swing Head**

Take care to ensure that the Swing Head is connected before the 2.855.0010 and 2.855.0020 instrument versions are set to work. Check the connection cable.

The connection socket (Mini DIN) for the Swing Head drive is located on the rear of the tower above the stirrer connector.



Figure 8 Connecting the Swing Head

If the Swing Head is not connected, connect it as follows:

1 Plugging in the cable

- Guide the connection cable of the Swing Head through the guide chain of the tower (*see chapter 3.9, page 27*).
- Plug the Mini DIN plug into the socket 'Swing Head'.

Tower configuration

Make sure that the correct axial distance is entered when configuring the tower. The axial distance depends on whether the 855 Robotic Titrosampler is used with or without Swing Head. The correct axial distance is shown in the following table.

Table 1Axial distance

Operation	Axial distance
with Swing Head	196 mm
without Swing Head	166 mm

3.5 Configuring the robotic arm

A variety of differently constructed robotic arms is available for enabling a wide range of applications. They differ from one another in their geometric sizes, e.g. swing radius or maximum permissible swing angle. The configuration data must be entered in the Sample Processor or in the control software **prior to the assembly of the robotic arm**. The data required is engraved on the underside of the robotic arm. Examples of the most common robotic arms are shown in the illustration below.



Figure 9 Robotic arms - standard model versions:

- **1 Transfer robotic arm (6.1462.030)** For sample transfer, left-swinging.
- **3 Titration robotic arm (6.1462.050)** With titration head, left-/right-swinging *).
- 5 Macro robotic arm (6.1462.070) With holder for a 6.1458.XXX titration head insert, right-swinging.
- 2 Transfer robotic arm (6.1462.040) For sample transfer, right-swinging.
- 4 Macro robotic arm (6.1462.060) With holder for a 6.1458.XXX titration head insert, left-swinging.

*) can be mounted in two ways



NOTICE

A detailed list of the available robotic arms, along with the necessary configuration data, can be found in *Chapter Robotic arms, page 51ff.*

The following figure illustrates the most important configuration data that needs to be set in the control software to ensure correct usage of a robotic arm (left-swinging, here).



Figure 10 Configuration data of the robotic arms

1 Swing axis

This runs through the middle of the Swing Head drive.

3 Source axis

This runs from the swing axis to the midpoint of the sample rack and marks the initial position of the robotic arm.

5 Max. swing angle

This stands for the swing range that the robotic arm can reach. The range runs from the source axis to the maximum possible robotic arm position.

2 Swing radius

This is determined by the length of the robotic arm. The radius runs from the axis of rotation to the midpoint of the tip of the robotic arm.

4 Swing offset

This determines the 0° position of the robotic arm.

Swing direction

Left-swinging (**swing direction +**) or right-swinging (**swing direction** –) model versions are available as different types of robotic arms. Left-swing-ing means swinging from the initial position (pointing towards the middle of the rack) outwards to the left.

In the case of a Sample Processor with two towers, a right-swinging robotic arm must be mounted on Tower 1, a left-swinging robotic arm on Tower 2. If the alignment is incorrect, the two robotic arms could possibly come into contact with one another, resulting in damage to the drives.

3.6 Mounting the robotic arm

Depending on the model, robotic arms can be mounted as either rightswinging or left-swinging. The position of the limitation screw of the robotic arm must be taken into account during assembly. The limitation screw must face the tower of the Sample Processor when mounting the robotic arm. The following illustration shows on the left the position of the limitation screw on a right-swinging robotic arm (**Position 1**) and on the right on a left-swinging robotic arm (**Position 2**).



Figure 11 Limitation screw on the robotic arm

For robotic arms which can be mounted in two different ways (e.g. 6.1462.050), the limitation screw can be fitted in accordance with the required assembly direction (see above).



NOTICE

The configuration data of a robotic arm must be configured in the control software **before it is mounted** (see chapter 3.5, page 19).

Mounting the robotic arm

The mounting of a robotic arm on the Swing Head is described here, taking as an example a 6.1462.070 robotic arm for titration and a 6.1458.040 titration head insert. Initialize the Sample Processor before mounting. After the initialization of the Sample Processor, the drive disk of the Swing Head is positioned as though the robotic arm were located in the outermost position.



Figure 12 Mounting the robotic arm Mount the robotic arm as follows:

- 1 Place the **6.1458.040 titration head insert** in the opening of the robotic arm. The marking on the edge must be closer to the Swing Head. Tighten the titration head insert with the provided screws.
- 2 Hold the robotic arm in such a way that the opening faces to the right. While doing so, rotate the robotic arm outwards as far as possible, i.e. towards the tower see above. Slip the robotic arm from below over the guide bolts of the drive disk of the Swing Head.



NOTICE

Take care not to twist the drive disk and thereby put strain on the drive.

3 Tighten the robotic arm to the Swing Head with the screws and washers supplied.

3.7 Robotic arms with beaker sensor

For safety reasons, the presence of a beaker on the sample rack of a Sample Processor can be detected. Some robotic arm model versions are therefore equipped with a beaker sensor.



Figure 13 Connecting a beaker sensor (for example 6.1462.150)

1 Plug of the connection cable

2 Connection socket on the 786 Swing Head

3 Beaker sensor

Contact sensor in accordance with the Piezo principle

A robotic arm with beaker sensor is mounted as described on *page 21*. The connection of the sensor cable must take place while the instrument is switched off.

The beaker sensor is automatically recognized when switching on the instrument.

Functioning of the beaker sensor

If the beaker sensor of the robotic arm is activated, then the lift of the Sample Processor will move automatically into its work position after a MOVE command. The presence of the sample vessel is checked by the robotic arm setting down on top of it. No separate LIFT command is required in such cases.



NOTICE

The work position of the lift must be configured in such a way that the robotic arm is in place on the sample vessel. The robotic arm must bend very slightly while doing so, so that the Piezo sensor will generate a signal.

3.8 Installing the rinsing and aspiration equipment

Various tubings are necessary for rinsing the electrode and the dosing tips as well as for aspirating the sample solution after the titration. First, mount the tubings on the distributor.

Mounting the rinsing and aspiration tubings



Install the tubings as follows:

Figure 14 Mounting the rinsing and aspiration tubings

1 Mounting the rinsing tubing

 Manually tighten the three 6.1805.060 FEP tubings (60 cm) in the M6 bore holes of the distributor. Place the tubings into the guide chain (see chapter 3.9, page 27). These are the feed lines for the spray nozzles.

2 Mounting the aspiration tubing

Manually tighten the 6.1805.510 FEP aspiration tubing (60 cm) in the M8 bore hole of the distributor.

3 Mounting the feed line for the rinsing liquid

Remove the union nut of the left-hand connector of the distributor and guide it over the end of a 6.1812.000 PTFE tubing.
 You may have to widen the tubing end in order to be able to better mount the tubing, see note below. Pull the end of the tubing over the connection nipple of the distributor and fasten in place with the union nut.

The tubing leads to the rinsing pump (**Pump 1**) and can be cut to the correct length.

NOTICE

The opening of the tubing may need to be widened with a pointed object (e.g. with a Phillips screwdriver).

A piece of sandpaper may be used to get a better grip on the tubing.

Do not widen the tubing end before having slid the union nut onto the tubing.

4 Mounting the outlet tubing

Remove the union nut of the right-hand connector of the distributor and guide it over the end of the 6.1812.000 PTFE tubing.
 Pull the end of the tubing over the connection nipple of the distributor and fasten in place with the union nut.
 The tubing leads to the aspiration pump (Pump 2) and can be cut

to the correct length.

Mounting the distributor

Proceed as follows to install the 6.1808.170 distributor on the guide chain:

1 Attaching the distributor

• Insert the bottom part of the distributor holder between two clips of the chain links. The tubing can remain connected.



• Insert the upper hook of the holder in the upper clip.



• Let the bottom part of the holder snap in the bottom clip.



2 Fixing the rinsing tubings

• Place the rinsing tubings into the guide chain (see chapter 3.9, page 27).

3.9 Placing tubing and cables in the guide chain

Tubing and cables can be placed in the guide chain.

The guide chain contains a firmly installed clip on each chain link.



CAUTION

When mounting tubing and cables, make sure that there is no traction on the drives while moving the lift or swinging the robotic arm. Traction on the drive can overload and damage the drive.

If a robotic arm is used, we recommend placing the tubing and cables only above the third chain link in the guide chain to prevent traction on the drives.

Do **not** or only partially place rigid tubings, as for example aspiration tubings made of PTFE, into the guide chain.



NOTICE

Make sure that tubing and cables do not kink.

Inserting and removing the tubing and cables requires no tools.

1 Placing tubing and cables

• Press one side of the clip downwards and place tubing or cables in the guide chain.



2 Removing tubing and cables

• Press one side of the clip downwards and remove tubing or cables from the guide chain.

3.10 Equipping the titration head

Mounting the rinsing and aspiration tubings

Example of a 6.1462.070 robotic arm with 6.1458.040 titration head.


Figure 15 Installing the rinsing tubings and the aspiration tip Proceed as follows:

1 Connecting the rinsing nozzles

• Connect the three rinsing tubings that are already connected to the distributor on the tower to the preinstalled rinsing nozzles on the titration head.

2 Inserting the aspiration tip

 Insert the 6.1543.170 aspiration tip into the left-side opening on the front of the titration head (only with 6.1458.040 titration head)

3 Connecting the aspiration tubing

• Connect the 6.1805.510 aspiration tubing already connected to the distributor with the aspiration tip.



Inserting the stirrer and the electrode, connecting the dosing tubings

Figure 16 Installing stirrer and electrode

The equipment of the titration head is completed as follows:

1 Inserting the rod stirrer

- Insert the rod stirrer (802 Stirrer) into the rear opening of the titration head (at the arrow).
- Insert the cable into the guide rail.

2 Mounting the stirring propeller

• Fasten the 6.1909.0XX stirring propeller to the rod stirrer from below.

3 Inserting the electrode

Insert the electrode (e.g., a 6.0229.100 Solvotrode) with a 6.1236.020 SGJ sleeve into the titration head.

4 Connecting the dosing tubings

Connect two dosing tubings (e.g., 6.1805.120) to the pre-mounted dosing tips on the titration head.

3.11 Connecting the tower stirrer

A DIN socket for connecting a rod stirrer (**802 Stirrer**) or a magnetic stirrer (**741 Stirrer**) is located on the rear side of the tower.



Figure 17 Rod stirrer 802 Stirrer



Figure 18 Magnetic stirrer 741 Stirrer

Take care to observe correct orientation of the contact pins when plugging in the stirrer connection cable. The rib on the outside of the plug must match the reference mark (on the left) on the socket.



Figure 19 Connecting the tower stirrer



NOTICE

If an MSB stirrer is connected to the **MSB1** socket, then the stirrer connector at tower 1 cannot be used, because both sockets are controlled internally via MSB1.

3.12 Connecting an external pump

The model version with built-in pump allows to connect an external pump at the rear of the tower. To aspirate clear aqueous media, the **823 Membrane Pump Unit** can be used. To aspirate organic media or aqueous media containing solids, the peristaltic pump **772 Pump Unit** must be used.

Connecting a pump



Figure 20 Connecting a pump

Connect the external pump as follows:

- Plug the threaded plug of the connection cable into the connection socket Ext. Pump 2 on the rear of the tower.
 Pay attention to the orientation of the three contact pins.
 - Tighten the knurled screw at the front end of the plug by hand in clockwise direction.



This secures the plug.

3.13 Connecting MSB devices

In order to connect MSB devices, e.g. stirrers or dosing devices, Metrohm instruments are equipped with up to a maximum of four connectors on what is referred to as the *Metrohm Serial Bus* (MSB). Various kinds of peripheral devices can be connected in sequence (in series, as a "Daisy Chain") at a single MSB connector (8-pin Mini DIN socket) and controlled simultaneously by the respective control instrument. In addition to the connection cable, stirrers and the Remote Box are each equipped with their own MSB socket for this purpose.

The following figure provides an overview of the instruments that can be connected to an MSB socket, along with a number of different cabling variations.



Dosimat / Dosino



The control instrument determines which peripheral devices are supported.



When connecting MSB devices together, the following must be observed:

- Only one device of the same type can be used at a single MSB connector at one time.
- Dosing devices of the 700 Dosino and 685 Dosimat plus type cannot be connected together with other MSB devices on a shared connector. These dosing devices must be connected separately.



CAUTION

Exit the control software before you plug in MSB devices. When it is switched on, the control instrument automatically recognizes which device is connected to which MSB connector. The operating unit or the control software enters the connected MSB devices into the system configuration (device manager).

MSB connections can be extended with the 6.2151.010 cable. The maximum connection length permitted is 15 m.

3.13.1 Connecting a dosing device

Three dosing devices can be connected to the instrument.

The types of dosing devices that are supported are:

- 800 Dosino
- 700 Dosino
- 805 Dosimat
- 685 Dosimat plus



WARNING

If a Dosino is connected to the 855 Robotic Titrosampler, then the connection cable must be equipped with a T.2400.102 ferrite core. The ferrite core diminishes any interference voltages and thus ensures compliance with the strict EMC standards in accordance with the applicable technical standards, see EU declaration of conformity of the instrument.

Proceed as follows:

1 Mounting the ferrite core

• Fasten a T.2400.102 ferrite core to the Dosino connection cable close to the plug.

2 Connecting a dosing device

- Exit the control software.
- Connect the connection cable of the dosing device to one of the sockets marked with **MSB** on the rear of the control device.
- Start the control software.



Figure 22 Connecting a dosing device

3.13.2 Connecting a stirrer or titration stand

You can use the following instruments:

These devices have a built-in magnetic stirrer (stirring "from below"):

- 801 Stirrer
- 803 Ti Stand

This device has no built-in magnetic stirrer (stirring "from above"):

• 804 Ti Stand with rod stirrer 802 Stirrer

Connect a stirrer or a titration stand as follows:

1 Connecting the stirrer or titration stand

- Exit the control software.
- Connect the connection cable of the magnetic stirrer or of the titration stand to one of the sockets marked with **MSB** on the rear of the control instrument.
- 804 Ti Stand only: Connect the rod stirrer to the stirrer connector (socket with stirrer symbol) of the titration stand.
- Start the control software.



Figure 23 Connecting an MSB stirrer







If an MSB stirrer is connected to the **MSB1** socket, then the stirrer connector at tower 1 cannot be used, because both sockets are controlled internally via MSB1. In case of a 2-tower model of a USB Sample Processor, this also applies for the MSB2 socket and the stirrer connector on tower 2.

3.13.3 Connecting a Remote Box

Instruments that are controlled via remote lines and/or that send control signals via remote lines can be connected via the 6.2148.010 Remote Box. In addition to Metrohm, other instrument manufacturers also use similar connectors that make it possible to connect different instruments together. These interfaces are also frequently given the designations "TTL Logic", "I/O Control" or "Relay Control" and they generally have a signal level of 5 volts.

Control signals are understood to be electrical line statuses or electrical pulses (> 200 ms) which display the operating status of an instrument or which trigger or report an event. Sequences on a variety of instruments can thus be coordinated in a single complex automation system. However, no exchange of data is possible.

Proceed as follows:

1 Connecting the Remote Box

- Exit the control software.
- Connect the Remote Box connection cable to one of the sockets marked with **MSB** on the rear of the control instrument.
- Start the control software.



Figure 25 Connecting the Remote Box

You can connect the following instruments to the remote connector, among others:

- 849 Level Control (fill level monitoring in a canister)
- 731 Relay Box (switch box for 230/110 volt alternating current sockets and low-voltage direct current outlets)
- 843 Pump Station (for complex sample preparations or for rinsing external titration vessels)

The Remote Box also has an MSB socket at which a further MSB device, e.g. a dosing device or a stirrer, can be connected.

You will find precise information concerning the pin assignment of the interface on the Remote Box in the *appendix*.

3.14 Connecting USB devices

Two USB connectors (type A sockets) are available for connecting instruments with an USB interface. The 855 Robotic Titrosampler then functions as USB hub (distributor). If you wish to connect more than two USB devices, you can also use an additional, commercially available USB hub.



NOTICE

When connecting an USB device, the control instrument automatically recognizes which device is connected. The control software enters a connected USB device automatically into the system configuration (device manager).

3.14.1 Connecting a barcode reader

A barcode reader is used as an aid for text and numerical input. You can connect a barcode reader with USB interface.

Connect the barcode reader as follows:

1 Connecting the cables

 Insert the USB plug (type A) of the barcode reader into one of the USB sockets on the rear of the instrument.



Figure 26 USB connectors

2 Configuring the barcode reader in the control software

• Configure the barcode reader in the configuration of the control software, as described in the online help of the software.

Setting the barcode reader

The barcode reader requires certain basic settings. You can find instructions on how to program the barcode reader in the manual for the barcode reader. Set the barcode reader into program mode and configure the following settings:

- Specify the keyboard layout for the desired country (USA, Germany, France, Spain, German-speaking Switzerland). This setting must match the setting in the control software.
 - Make sure that the Ctrl characters (ASCII 00 to 31) can be sent.
 - Set the "Preamble" (introduction) or "Prefix Code" for the ASCII character 02 (STX or Ctrl B) to be sent as the first character.
 - Set the "Postamble", "Record Suffix" or "Postfix Code" for the ASCII character 04 (EOT or Ctrl D) to be sent as the last character.
 - Close programming mode.

3.15 Mounting the base plate

If it is needed by an application that a determination is not carried out on the sample rack but in an external measuring cell, a stand plate can be mounted. It can be placed on the left or on the right of a tower of the 855 Robotic Titrosampler. Any accessories parts can be placed on the support rod of the stand plate.

Mounting the base plate



Figure 27 Mounting the base plate

Proceed as follows:

- 1 Insert the enclosed countersunk screw from below into the opening on the very back of the base plate.
- 2 Place the washer (flat side upwards) over the countersunk screw and tighten the support rod with the countersunk screw. The necessary hexagon key is enclosed with the 855 Robotic Titrosampler.
- **3** Hang the base plate to the assembly rail of the turntable with the holding clamps. Slide the whole stand plate as near to the tower as possible. Fix the holding clamps with a hexagon key to the assembly rail.

i

NOTICE

Before the stand plate is fixed to the tower, e.g. a magnetic stirrer or a measuring cell can be mounted on the support rod.

- **4** Guide the bracing from above over the support rod.
- **5** Loosen the screw for fastening the bracing on the rear panel of the tower (see figure) and fix the bracing with the enclosed hexagon screw.
- 6 Fix the bracing to the support rod with a hexagon key, see figure.

3.16 Mounting the drip pan

Serious damage to the instrument or a danger to the user can occur if chemicals or liquid samples are spilled.

The use of a suitable **drip pan** (6.2711.0X0) is recommended in order to avoid such incidents.

Mounting the drip pan



Figure 28 Installing the drip pan Install the drip pan as follows:

- **1** Fasten the enclosed tubing to the drainage nipple on the drip pan and then guide the free end of the tubing into a waste container.
- **2** Place the drip pan on the assembly rail of the turntable as shown in the figure.

3.17 Attaching a sample rack



Figure 29 Attaching rack

The easiest way to put a sample rack into position is when the turntable is in the starting position.

If the instrument is switched off, then the turntable can be rotated manually into position. Both of the turntable guide bolts must be positioned so that they are pointing towards the tower.

Put the rack into place as follows:

1 Carefully center the rack on the turntable. The guide bolts on the turntable must engage with the openings in the bottom of the rack.



NOTICE

Hold the rack in such a way that the printed Metrohm logo is legible horizontally.

2 Screw the fixing screw in the handle tight by rotating it clockwise.



3 Carry out the **[Rack reset]** or **[Initialize rack]** function in the manual control of the control software.

The rack is moved into starting position. The magnet code of the rack is read out by the instrument during this process. The white arrow in *Figure 29* Figure indicates the position of the magnet holder. The six-

digit magnet code is used to identify the rack type. The sample positions and any special positions on the rack are defined along with the rack type.

3.18 Mounting the safety shield



WARNING

It is imperative that the safety shield be installed before the first time the 855 Robotic Titrosampler is used. The device may not be operated without a safety shield.



Figure 30 Mounting the safety shield Proceed as follows.

- **1** Loosen the knurled screws on both sides of the tower.
- **2** Move the **safety shield** into position, starting from the top. Observe the corresponding illustration at the top.
- **3** Fix the safety shield in place with the knurled screws.

i

NOTICE

You can adjust the vertical position of the safety shield at any time by loosening the screws.

The optimum vertical position meets the following criteria:

- The rack can rotate freely.
- Nobody can access the working area of the lift while the instrument is in operation.

4 Operation and maintenance

4.1 General notes

The 855 Robotic Titrosampler requires appropriate care. Excess contamination of the instrument may result in malfunctions and a reduction in the service life of the sturdy mechanical and electronic components.

Severe contamination can also have an influence on the measured results. Regular cleaning of exposed parts can prevent this to a large extent.

Spilled chemicals and solvents must be removed immediately. In particular, the power plug should be protected from contamination.

4.2 Care

- Check all tubing connections regularly for leaks.
- Rinse the tubing connections from time to time. The tubing must be replaced after prolonged usage.

5 Troubleshooting

5.1 Sample Processor

Problem	Cause	Remedy
The instrument is not recognized.	Sample Processor – No USB connection.	 Correctly plug in the USB connecting cable on both ends. Restart the control software or switch the Touch Control off and on again.
	Sample Processor – Power supply of the instrument is missing.	 Plug in the power cord on the instrument. Restart the control software or switch the Touch Control off and on again.

5.2 Robotic arm

Problem	Cause	Remedy
The robotic arm moves all the way outward and buz- zes.	Sample Processor – The Swing Head is not correctly configured.	In the control software under "Configuration" (or under "Device manager" for Touch Con- trol), enter the correct value for the Swing offset .
	Sample Processor – Robotic arm is mounted incorrectly.	Disconnect the power plug and dismount the robotic arm. Check the configuration of the robotic arm and mount it correctly if necessary (left-swinging ⇔ right-swinging).
The Swing Head either misses the rack positions totally or is inaccurate	Sample Processor – The Swing Head is not correctly configured.	Enter the correct values for the Swing radius , Swing offset etc. in the control software under "Configuration" (or under "Device man- ager" for Touch Control).
	Sample Processor – The axial distance is not cor- rectly configured.	Enter the correct value for the Axial distance in the control software under "Configuration" (or under "Device manager" for Touch Con- trol).
	Sample Processor – The wrong rack table is being used.	Initialize the rack using the function Initialize rack in the "Manual control".

Problem	Cause	Remedy	
	Swing Head – The Swing Head drive is defective.	Contact Metrohm Service.	

5.3 Pump

Problem	Cause	Remedy
The pump is leaking.	Sample Processor – A tub- ing connection is leaking.	Check the tubing connections especially between the distributor and the pump and seal tightly.
	Canister – There is too much pressure on the pump valve.	 Make sure that the canisters are not placed on a higher level than the pump. Check the fill level of the canisters. Make sure that the canisters are not com- pletely sealed (overpressure).

6 Appendix

6.1 Beaker sensor

Every tower of a Sample Processor is equipped with a beaker sensor detecting the availability of a sample vessel in front of the tower. An infrared sensor identifies devices of various materials if they are located in a correct position in front of the tower. In the rack configuration of the control instrument or the control software, **Tower** must be selected for the setting 'Beaker sensor'. This beaker test is carried out whenever a rack position is approached in a method run.



Figure 31 Beaker sensor on the tower

The beaker sensor on the tower can only be used with single-row sample racks.

6.2 Rinsing nozzles

Using rinsing nozzles is very effective in order to rinse sample vessels (with sensors and buret tips) efficiently. Rinsing nozzles are available in two model versions:

6.2740.020 spray nozzle

For the fine-spraying of the rinsing solution. The nozzle has a small ball at the opening. The distribution (but also the backpressure) of the rinsing liquid is clearly higher than the one of a rinsing nozzle.

• **6.2740.030 rinsing nozzle** (optionally available) The rinsing liquid is applied as a fine jet for optimal removal of layers on electrodes and on titration accessories.



Figure 32 Spray nozzles - Mode of operation

The height of the nozzles can be adjusted in the titration head in order to reach an optimal rinsing effect.

6.3 Remote interface

The 6.2148.010 Remote Box allows devices to be controlled which cannot be connected directly to the MSB interface of the Titrosampler.



Figure 33 Connectors of the Remote Box

1 Cables

For connecting to an MSB connector of the Titrosampler.

3 Remote connector

For connecting instruments with a remote interface.

2 MSB connector

Metrohm Serial Bus. For connecting external dosing devices or stirrers.

6.3.1 Pin assignment of the remote interfaces



Figure 34 Pin assignment of remote socket and remote plug

The above figure of the pin assignment of a Metrohm remote interface does not only apply for the Remote Box, but also for all Metrohm devices with 25-pin D-Sub remote connectors.

Inputs



approx. 50 k Ω Pull-up $t_p > 20 \text{ ms}$ active = low, inactive = high

The input lines can be scanned with the **SCAN** command.

Outputs



The output lines can be set with the **CTRL** command.

Table 2Inputs and outputs of the remote interface

Pin no.	Assignment	Pin no.
21	Output 0	5
9	Output 1	18
22	Output 2	4
10	Output 3	17
23	Output 4	3
11	Output 5	16
	21 9 22 10 23	21 Output 0 9 Output 1 22 Output 2 10 Output 3 23 Output 4

1
2
6
7
8
13
19
20
-

Stirring rate 6.4

The stirring rate can be adjusted in steps from -15 to +15.

The approximate rotational speed for the internal magnetic stirrer (depends on the product version) can be calculated with the following formula:

Rotational speed/min
$$(r/min) = 125 \cdot Stirring$$
 rate

Example:

Configured stirring rate: 8

Rotational speed in revolutions per minutes = $125 \cdot 8 = 1,000$



The information on the separately connectable 802 propeller stirrer can be found in the "802 Stirrer" manual.

6.5 Robotic arms

6.5.1 Robotic arms for titration

Table 3Configuration data of the titration robotic arms

Туре	6.1462.050	6.1462.060	6.1462.070	6.1462.260
	6.1462.050	6.1462.060	6.1462.070	6.1462.260
Swing direction	+/	+	_	+
Swing offset	0°	-8°	-8°	-8°
Max. swing angle	84°	73°	73°	105°
Swing radius	110 mm	127 mm	127 mm	110 mm

6.1462.050	Robotic arm with titration head, left or right-swinging
	For titration in 75 mL sample vessels and larger.
	The arm can be equipped with two microelectrodes, one propeller stirrer and three spray nozzles. Two buret tips with anti-diffusion valve and one aspiration tip with connections for M6 tubing are already retracted into the arm.
	Material: PP
6.1462.060	Robotic arm with holder for a titration head, left-swinging
	The arm can be modified to create the desired titration robotic arm by means of the installation of a titration head 6.1458.xxx.
	Material: PP
6.1462.070	Robotic arm with holder for a titration head, right-swinging
	The arm can be modified to create the desired titration robotic arm by means of the installation of a titration head 6.1458.xxx.
	Material: PP

6.1462.260 Robotic arm with holder for a titration head, left-swinging, external

The arm can be modified to create the desired titration robotic arm by means of the installation of a titration head 6.1458.xxx. The cutout permits movement to external positions near the rack, e.g. an external rinsing station.

Material: PVC

Туре	6.1462.150	6.1462.160	6.1462.170

Table 4 Configuration data of the robotic arms with beaker sen	Table 4	Configuration	data	of the	robotic	arms	with	beaker sensc	r
--	---------	---------------	------	--------	---------	------	------	--------------	---

Туре	6.1462.150	6.1462.160	6.1462.170
	6.1462.150	6.1462.160	6.1462.170
Swing direction	+/	+	-
Swing offset	0°	-8°	-8°
Max. swing angle	84°	73°	73°
Swing radius	110 mm	127 mm	127 mm

6.1462.150 Robotic arm with titration head and beaker sensor, left or rightswinging

For titration in 75 mL sample vessels and larger.

The arm can be equipped with two microelectrodes, one propeller stirrer and three spray nozzles. Two buret tips with anti-diffusion valve and one aspiration tip with connections for M6 tubing are already retracted into the arm.

Material: PP

6.1462.160 Robotic arm with holder for a titration head and beaker sensor, left-swinging

The arm can be modified to create the desired titration robotic arm by means of the installation of a titration head 6.1458.xxx.

Material: PP

6.1462.170 Robotic arm with holder for a titration head and beaker sensor, right-swinging

The arm can be modified to create the desired titration robotic arm by means of the installation of a titration head 6.1458.xxx.

Material: PP

6.5.2 **Robotic** arms for sample preparation

Table 5	Configuration	data of the robotic	arms for sample p	preparation	

Туре	6.1462.030	6.1462.040	6.1462.090	6.1462.240
	6.1462.030	6.1462.040	6.1462.090	6.1462.240
Swing direction	+	-	-	-
Swing offset	8°	8°	8°	-8.6°
Max. swing angle	117°	117°	117°	122°
Swing radius	112 mm	112 mm	112 mm	149.8 mm

Robotic arm with transfer head, left-swinging
Robotic arm for fully automated pipetting or dilution of liquid samples with Sample Processor Systems.
Material: PP
Robotic arm with transfer head, right-swinging
Robotic arm for fully automated pipetting or dilution of liquid samples with Sample Processor Systems.
Material: PP
Robotic arm with Luer lock adapter, right-swinging
For the connection of hollow needles with Luer lock connection. Suitable for the transfer of samples from sealed vials with septum seal.
Material: PP
Robotic arm with transfer head, bent, right-swinging
The transfer head can, when equipped with 6.1808.220 adapter, be used as a holder for various tools with Luer connection on multirow racks.

Material: PP

6.5.3 **Robotic arms for special applications**

Table 6	Configuration	data of the special robotic arms
---------	---------------	----------------------------------

Туре	6.1462.250	6.1462.080
	6.1462.250	6.1462.080
Swing direction	+	+
Swing offset	0°	8°
Max. swing angle	115.5°	117°
Swing radius	110 mm	112 mm

6.1462.250 Robotic arm as holder for a Polytron, left-swinging

The robotic arm makes it possible to use the Polytron for sample preparation on multirow sample racks. It contains one retracted buret tip for adding solvents and three spray nozzles for cleaning.

Material: PP

6.1462.080 Robotic arm DIS-COVER, left-swinging

Robotic arm for placing and removing sample vessel covers (75 and 250 mL) covers 6.2037.050 und 6.2037.060) on the sample rack of a Robotic Sample Processor.

Material: PP

7 Technical specifications

7.1 Measuring interface

The 855 Robotic Titrosampler has one galvanically isolated measuring interface.

The measuring cycle is 100 ms for all measuring modes.

7.1.1 Potentiometry

One high-ohm measuring input (**Ind.**) for pH, redox or ISE electrodes and one reference measuring input (**Ref.**) for a separate reference electrode.

Input resistance	$> 1 \cdot 10^{12} \Omega$
Offset current	$< 1 \cdot 10^{-12}$ A (under reference conditions)
Measuring mode	
Measuring range	–20 - +20 pH
Resolution	0.001 pH 0.1 mV
Measuring accu- racy	±0.003 pH ±0.2 mV (±1 digit, without sensor error, under reference conditions)

7.1.2 Polarizer

	A measuring input (Pol.) for a polarizable electrode.
Measuring mode Ipol	Determination with adjustable polarization current
Polarization cur- rent	$-125.0\ldots+125.0~\mu A$ (in steps of 2.5 μA) -125 - $-121~\mu A$ / $+121$ - $+125~\mu A$: non-guaranteed values, dependent on reference voltage +2.5 V
Measuring range	-1200 - +1200 mV
Resolution	0.1 mV
Measuring	±0.2 mV
accuracy	(±1 digit, without sensor error, under reference conditions)
Measuring mode Upol	Determination with adjustable polarization voltage
Polarization voltage	-1250 - +1250 mV (in steps of 25 mV)

	–1250 - –1210 mV / +1210 - +1250 mV: non-guaranteed values, dependent on reference voltage +2.5 V
Measuring	–120 - +120 μA
range	
Resolution	0.1 μΑ

7.1.3 Temperature

One measuring input (**Temp.**) for a Pt1000 or NTC temperature sensor. Automatic temperature compensation for NTC sensors; R (25 °C) and B value can be configured.

Measuring range	
Pt1000	–150 - +250 °C
NTC	–5 - +250 °C
	(R (25 °C) = 30,000 Ω and B (25/50) = 4,100 K)
Resolution	
Pt1000	0.1 °C
NTC	0.1 °C
Measuring accu-	
racy	
Pt1000	±0.2 °C
	(Applies for measuring range $-20 - +150 \text{ °C}$; $\pm 1 \text{ digit}$; without sensor error, under reference conditions)
NTC	±0.6 °C (Applies for measuring range +10 - +40 °C; ±1 digit; without sensor error, under reference conditions)

7.2 Lift and turntable

Stroke path	235 mm
Maximum lift load	approx. 30 N / 3 kg
Lift rate	adjustable, 5 - 25 mm/s
Turntable speed	adjustable, 3 - 20 angular degrees/s

7.3 Membrane pump(s) with valve

Flow rate	> 450 mL/min
	Pressure head 2 m

7.4 786 Swing Head

Maximum load	approx. 15 N
Swing rate	10 - 55 angular degrees/s
Beaker sensor	Socket with M8 thread
connector	

7.5 Interfaces and connectors

Connector "Con- troller"	USB upstream port (9-pin Mini DIN socket) for connecting a computer for controlling the instrument.
MSB connectors MSB1 - MSB3	Three 9-pin Mini DIN sockets for connecting dosing devices (Dosino/ Dosimat), stirrers, etc.
USB connectors 1/2	Two USB downstream ports (type A sockets), each 500 mA, for con- necting Metrohm instruments or USB peripheral devices of other man- ufacturers.
Stirrer connector	DIN socket
Stirring rate	722/802 rod stirrer: 180 - 3000 rpm 741 magnetic stirrer: 180 - 2600 rpm adjustable in 15 steps each in both directions of rotation
Pump connector	Two sockets with M8 thread for 772 Pump Unit or 823 Membrane Pump Unit U = $16 \pm 1 \text{ V}$, I = $\leq 0.8 \text{ A}$
Swing Head con- nector	9-pin Mini DIN socket

7.6 **Power connection**

Voltage	100 - 240 V
Frequency	50 - 60 Hz
Power consump- tion	115 W
Fuse	2.0 ATH

7.7 Ambient temperature

Nominal function	5 - 45 °C
range	Humidity < 80%
Storage	–20 - 60 °C
Transport	–40 - 60 °C

7.8 Reference conditions

Ambient tempera-	25 °C (±3 °C)
ture	
Relative humidity	$\leq 60\%$

7.9 Dimensions

Width	0.28 m
Height	0.73 m
Depth	0.53 m
Weight (without accessories)	1.855.0010: 15.5 kg 1.855.0020: 16.4 kg 1.855.0030: 15.5 kg 1.855.0040: 16.4 kg
Material	
Housing	Metal housing, surface-treated

8 Accessories

Up-to-date information on the scope of delivery and optional accessories for your product can be found on the Internet. You can download this information using the article number as follows:

Downloading the accessories list

1	Enter https://www.metrohm.com/ into your Internet browser.
2	Enter the article number (e.g. 855) into the search field. The search result is displayed.
3	Click on the product.
	Detailed information regarding the product is shown on various tabs.
4	On the Included parts tab, click on Download the PDF . The PDF file with the accessories data is created.
	NOTICE
On	ce you have received your new product, we recommend download-

Once you have received your new product, we recommend downloading the accessories list from the Internet, printing it out and keeping it together with the manual for reference purposes.

Index

Index

6.2151.000 controller cable	16
685 Dosimat plus	34
700 Dosino	34
741 Stirrer	31
800 Dosino	34
801 Stirrer	35
802 Stirrer	31
803 Ti Stand	35
804 Ti Stand	35
805 Dosimat	34

Α

Aspiration equipment	24
Axial distance	18

В

Back panel	12
Barcode reader	
Connect	37
Beaker sensor 23,	52

С

Computer
Connect 16
Connect
Barcode reader 37
Computer 16
Dosing device
MSB devices 33
Power grid 15
Remote Box 36
Stirrer 35
Titration stand
USB devices 37
Connector
MSB 4
USB 4
Connectors 12
Controller
Connector 12

D

Device software	
Update	2
DIS-COVER	. 54
Distributor	. 26
Dosing device	
Connect	. 34
Drip pan	40
Drive disk	. 22

Driver software	
Install	16

<u>E</u>
Electrode
Mount 30
Electrostatic charge 7

F Fe

errite core	
Mount	34

G

u .	
Guide bolts	22
Guide chain	27

I

Install	
Pump	32
Installation	
Driver software	16

L

Left-swinging	19
Limitation screw	21
Luer lock adapter	53

м

Magnet code 14
Magnetic stirrer
Mounting 31
Maintenance 44
Max. swing angle 20
Metrohm Serial Bus MSB, see also
"MSB"
Model versions 3
Mount
Aspiration equipment 24
Drip pan 40
Rinsing equipment
MSB
Connect devices
Connector 12
MSB connector 4

Ρ

Pins	49
Polytron	54
Power connection	12, 15, 16
Propeller stirrer	31

Pump	
Install	32

R

R
Rack 13
Rack code 14
Remote
Input 49
Interface 49
Output 49
Remote Box
Connect 36
Pin assignment 49
Right-swinging 19
Rinsing equipment 24
Robotic arm
Beaker sensor 23, 52
Bent 53
Configuration data 20, 51, 53,
54
Configuring 19
Left-swinging 21
Luer lock adapter 53
Model versions 19, 51
Mount 21
Right-swinging 21
Titration head 51
Transfer head 53
Robotic arms
Model versions 53
Rod stirrer 31

S

-
Safety instructions 6
Sample preparation 53
Sample rack 13
Serial number 12
Service 6
Stirrer
Connect 35
Connecting 31
Mount
Stirring rate 50
Supply voltage7
Swing axis 20
Swing direction 21
Swing Head
Connect 18
Swing offset 20

Swing radius 20	C
т	
Titration head 5	1
Equip 28	3
Titration stand	
Connect 3!	5

Tower configuration	18
Tower stirrer	31
Transfer head	53
U	
Update	

Device software 2

####