

BENCH-TOP SHAKER

Orbitron

A true workhorse

The extraordinarily stable, splash-proof Orbitron is suitable both for daily operations in the lab and for use in climatecontrolled rooms.

Up to 31 kg loading capacity

The Orbitron is the powerhouse of benchtop shakers and can be loaded with up to 31 kg of materials. Select from two different tray sizes. The clever quick-lock system lets you change trays in no time.

Lab or climate-controlled room? Both!

The Orbitron is a splash-proof shaker that is exceptionally well-suited both for daily operations in the lab and for use in climate-controlled areas.

Penta-Drive

The patented, high-performance Penta-Drive unit is at the heart of many INFORS HT shakers. The highlight of the drive is its ingenious dynamic mass-balancing system, which ensures whisper-quiet, vibration-free operation.



Shake flask sizes

25 mL - 5000 mL



Suitable for

Bacteria, algae and cell cultures



High-performance Penta-Drive



Compatibility

With the INFORS HT M- and N-trays, clamps, test tube holders, microtitre plates and Sticky Stuff.



Features & Accessories

Housing

• Rounded edges make it easy to clean

Drive

- Patented Penta drive with dynamic balancing
- Quiet and vibration-free
- Precise shaking speed control
- Watertight

Operating unit

- Splash-proof display (IP43)
- Displays shaking speed and time
- Integrated timer

Accessories

- Autoclavable trays made of corrosion resistant, anodized aluminium
- Sticky Stuff adhesive mats for a reliable fixation
- Stainless steel clamps for Erlenmeyer and Fernbach flasks
- Holders for test tubes or Falcon tubes

Technical data

Dimensions (W x D x H)	640 mm x 600 mm x 150 mm	
Number of batches	197	
Volume	21	
Maximum load	31 kg	
Rotation speed	20–550 min ⁻¹	
Shaking throw	25 mm	
Temperature	4 °C bis 65 °C	
Standard parameters	Rotation speed, timer	
Power supply	115/230 V ± 10 %, 50-60 Hz	

Tray capacity

Shake flask	M tray	N tray
size	(850 mm x 470 mm)	(480 mm x 420 mm)
25 mL	197	105
50 mL	131	69
100 mL	91	46
250 mL	48	25
500 mL	31	15
1000 mL	19	10
2000 mL	13	6
3000 mL	9	
4000 mL	8	
5000 mL	6	



«A simple shaker that shines for its reliable continuous operation, even with high loads.»

Damian von Blarer – Area Sales Manager



Contact us for a personal consultation.

www.infors-ht.com/en/orbitron

10.20 May be subject to technical amendments.

